

Kalmely.

THE CELLULAR RESET

The Cellular Reset Protocol

Your 3-minute nightly guide to calmer, clearer breathing — and evenings that finally feel simple again.

A GIFT FROM KALMEY · NOT FOR RESALE

This is the guide we wish someone had handed us

If a blocked-up nose has been quietly running your evenings — the ritual before bed, the little bottle within reach, the 3 a.m. wake-ups — you're in the right place.

Your Kalmely device uses gentle 660-nanometer red light, a spectrum studied in photobiomodulation research for decades. But the device is only half of it. **What changes your nights is the routine built around it** — and that's what this guide gives you.

Inside you'll find the science in plain English, the exact 3-minute nightly protocol, a 14-night calendar of what to expect, and the small evening habits that multiply your results.

HOW TO USE THIS GUIDE

Read Parts 1 and 2 tonight — fifteen minutes, tops. Start the protocol this evening. Keep the 14-night calendar (page 14) somewhere you'll see it. That's it.

What's inside

- 1 Why your nose feels worse at night — the swollen lining, explained

- 2 The quick-fix trap — what nobody tells you about decongestant sprays

- 3 The Cellular Reset — how 660nm red light works, in plain English

- 4 Your device — setup, fit, charging and hygiene

- 5 The 3-Minute Nightly Reset — the exact protocol

- 6 The 14-night calendar — what to expect, night by night

- 7 The Wind-Down Stack — five evening habits that multiply results

- 8 Retrain your breathing — two gentle exercises

- 9 Simplifying your routine — if you currently rely on sprays

- 10 Travel, care & troubleshooting

- 11 Safety notes & when to see a professional

PART ONE

Understand what's happening

Ten minutes of plain-English science — because a routine you understand is a routine you keep.

1

PART ONE · UNDERSTANDING

Why your nose feels worse at night

Here's the part almost nobody explains: the "blocked" feeling usually isn't mucus. It's the **lining** of your nose — soft, blood-rich tissue that swells and shrinks all day long.

That lining is one of the most reactive tissues in your body. Dry indoor air, dust, a glass of wine, hormonal shifts, even the natural rhythm of your own circulation — all of it can leave the lining puffier than it should be. When it swells, the airway narrows, and every breath feels like work.

Three reasons nights are the worst

- **Gravity changes.** When you lie down, more blood flows to the vessels in the nasal lining — so tissue that was only slightly puffy at dinner can feel fully sealed by midnight.
- **The room dries out.** Heating and air-conditioning strip moisture from the air, and a dry lining is an irritated lining.
- **The nasal cycle.** Your nostrils naturally take turns doing most of the breathing, alternating every few hours. On a calm lining you never notice. On a swollen one, each "off-duty" side feels completely shut.

THE KEY IDEA

If the swollen lining is the problem, then the goal isn't to force it open for a few hours. It's to help the tissue itself settle down — night after night.

2

PART ONE · UNDERSTANDING

The quick-fix trap

If you've ever felt like you *need* that little bottle on the nightstand — this chapter is for you, and there's nothing to be embarrassed about.

Decongestant sprays work by clamping the blood vessels in the nasal lining shut. The relief is real and fast. But when the effect wears off, the vessels often swing back — sometimes puffier than before. Doctors call it **rebound congestion** (rhinitis medicamentosa). You may know it as "it only works for a few hours now."

So the dose creeps up. The bottle multiplies — one in the purse, one in the car, one by the bed. Not because you did anything wrong, but because that's how the mechanism works. The spray's own label says "do not use for more than 3 days" for exactly this reason.

"It's not the congestion that wears you down. It's needing the bottle — and planning your whole evening around it."

The takeaway: anything that only forces the vessels shut is renting relief by the hour. A calmer lining doesn't need to be forced. In chapter 9 we'll cover how many people gradually simplify their spray routine — gently, and with their doctor or pharmacist in the loop.

3

PART ONE · UNDERSTANDING

The Cellular Reset: how 660nm light works

Your device doesn't push, spray or medicate anything. It delivers one very specific thing: **red light at 660 nanometers**, directly where the swollen lining lives.

Plain English, three steps

1

The light reaches the cells of the lining

660nm red light is absorbed by a molecule called cytochrome c oxidase, found in the mitochondria — the tiny power plants inside each cell.

2

The power plants get a nudge

With that light absorbed, mitochondria produce cellular energy (ATP) more efficiently — like a hard reset for cells that have been running hot and irritated.

3

The lining can settle

Research on this spectrum shows it helps modulate the signals associated with an irritated, swollen state — supporting the tissue's own return to calm, rather than forcing vessels shut.

WHERE THE SCIENCE COMES FROM

A double-blind trial published in *Annals of Allergy, Asthma & Immunology* (Neuman & Finkelstein, 1997) tested pure 660nm intranasal light: 72% of participants reported meaningful improvement, versus 24% with the placebo device. Individual results vary — but this is why we use this exact wavelength, and not a random red LED.

Why 3 minutes a night beats 30 minutes once

Photobiomodulation — the science of light and cells — follows a rule that surprises people: **small, regular doses outperform big, occasional ones.**

Cells respond to light the way muscles respond to exercise. One heroic gym session changes nothing; twenty short ones change everything. The lining of your nose renews and regulates itself on a cycle of days and weeks — so the win comes from showing up nightly, not from doubling the dose.

That's why the protocol is deliberately tiny. Three minutes, once a night, ideally at the same point in your evening. Small enough that you'll actually do it on the tired nights — and the tired nights are the ones that count.

SET YOUR EXPECTATION HERE

This is not a "feel it in 60 seconds" product — anything that works that fast is forcing the tissue, and you already know how that story ends. Most people describe the first differences as nights that feel *calmer* and mornings that feel less puffy, building over the first couple of weeks. Give it the full 14 nights before you judge it.

PART TWO

The Protocol

Your device, your 3-minute routine, and the 14-night calendar.

4

PART TWO · THE PROTOCOL

Your device, in two minutes

In the box

- The Kalmely device — the small oval unit with the control button.
- Two soft silicone light tips — these glow red and sit gently just inside each nostril.
- USB charging cable.

Before your first session

1

Charge it fully

Plug in the USB cable until the indicator shows a full charge. A full battery covers many sessions — you won't charge often.

2

Fit the tips

Seat each silicone tip firmly onto the device. They should feel snug, never loose.

3

Clean before first use

Wipe the tips with a little warm water and mild soap (or an alcohol wipe), and let them dry completely. Make this a weekly habit — hygiene first, always.

GOOD TO KNOW

The light should always look **red and steady**. Red light at this intensity is gentle, non-UV and doesn't heat the tissue — most people feel nothing at all, or the faintest sense of warmth.

5

PART TWO · THE PROTOCOL

The 3-Minute Nightly Reset

Do this once every evening, 30–60 minutes before bed. Same time, same place — attach it to something you already do (after brushing your teeth is the classic).

1

Settle somewhere comfortable

Sofa, bed, reading chair. You can read, watch TV or just breathe — this is a hands-free routine.

2

Place the tips

Gently rest a silicone tip just inside each nostril. Comfortable, never deep. If it feels intrusive, it's in too far.

3

Press the button

The tips glow red and the session begins. Breathe slowly through your nose if you can; through the mouth is fine too.

4

Let it finish

The session ends on its own after ~3 minutes. Remove the tips, set the device on the nightstand, done.

THE FIRST THREE NIGHTS

If you're sensitive, start with 1–2 minutes and work up to the full session. New sensations fade after a night or two — after that, most people say they forget it's on.

The three rules that make it work

Rule 1 — Never skip two nights in a row

Miss one night? Fine, life happens. But two missed nights become a lost week faster than you think. If you're exhausted, do a shortened one-minute session — the habit matters more than the dose.

Rule 2 — Anchor it, don't schedule it

Routines tied to a clock die in a week. Routines tied to another habit last years. "After I brush my teeth, I do my reset" — that phrasing, literally, is the whole trick.

Rule 3 — Track nights, not feelings

Day-to-day sensations bounce around — allergen exposure, dry air and sleep quality all add noise. Judge the trend at night 7 and night 14, not at 2 a.m. on night 3. Use the calendar on the next page.

"The people who write to us glowing at week three all say the same thing: they stopped overthinking it and just did their three minutes."

6

PART TWO · THE PROTOCOL

The 14-night calendar

Nights 1–7: building the habit.

Nights 1–2 · First contact

The novelty phase. You might feel gentle warmth, or nothing at all — both are normal. Your only job: do the session and put the device where you'll see it tomorrow.

Nights 3–4 · The honest stretch

Novelty fades and this is where routines usually die. Use the anchor habit. Some people notice evenings starting to feel a touch calmer here; many don't yet.

Nights 5–7 · First check-in

End of week one. Ask one question only: do the nights feel any calmer than last week? Don't grade individual nights — look at the week as a whole.

A TIP THAT HELPS

Jot one word per morning on your phone — "rough", "okay", "calm".
Fourteen words take no effort, and the pattern they reveal at the end is worth more than memory.

Nights 8–14: letting it build

Nights 8–10 · The quiet middle

The routine now runs on autopilot — about 20 seconds of actual effort. This is where consistent light exposure does its patient work on the lining. Keep the wind-down habits from Part Three going.

Nights 11–13 · Watch the mornings

For many people the first clear signal isn't at night — it's waking up with the nose feeling less puffy, or noticing the "off-duty" nostril doesn't feel completely sealed anymore.

Night 14 · The real verdict

Compare your fourteen morning words against week zero. Calmer trend? Keep going — most long-term users say weeks 3–6 are where it settled in for good. Flat trend? Check the troubleshooting list in chapter 10 before you conclude anything.

AFTER NIGHT 14

Keep the nightly session — this is a maintenance routine, not a course you finish. Three minutes is cheap insurance for calm nights.

PART THREE

Multiply your results

The evening habits, the breathing drills, and how to simplify the rest of your routine.

The Wind-Down Stack

Five small evening habits. None replaces your reset — but each one removes an irritant that works against a calm lining. Stack as many as you comfortably can.

1 • Humidify the bedroom

Dry air is the lining's worst roommate. A basic humidifier (aim for 40–50% humidity) or even a bowl of water near the radiator makes winter nights dramatically kinder to your nose.

2 • Steam before the session

A warm shower before your reset — or three slow minutes breathing over a warm mug — leaves the lining moist and comfortable. Lovely combination, especially in dry seasons.

3 • Mind the last glass

Alcohol dilates the blood vessels in the nasal lining — that's why wine can leave you stuffier within the hour. You don't have to give it up; just notice the pattern and keep evenings light when calm nights matter most.

The Wind-Down Stack, continued

4 · Elevate your head a little

Remember the gravity effect from chapter 1? An extra pillow — or raising the head of the bed a few centimeters — reduces the nighttime blood-pooling that swells the lining while you sleep. Side sleeping helps many people too.

5 · Rinse on high-exposure days

Pollen day? Dusty work? A simple saline rinse (plain salt water, from any pharmacy) before your evening reset physically washes irritants off the lining. Saline is not a decongestant spray — it's just clean water for your nose, and it pairs beautifully with the light.

THE ORDER THAT WORKS

Warm shower → saline rinse if needed → **3-minute reset** → bed with the humidifier running. Fifteen unhurried minutes, and every step is drug-free.

8

PART THREE · MULTIPLY

Retrain your breathing

After months of a blocked nose, many people quietly become mouth-breathers — and the nose, like anything unused, gets worse at its job. These two drills gently reverse that. Do them during your 3-minute session: perfect timing, zero extra effort.

Drill 1 • The slow six

Breathe in through the nose for a slow count of four, out through the nose for a count of six. Repeat for one minute. The extended exhale settles the nervous system — most people feel their shoulders drop by the third breath.

Drill 2 • The gentle hold

Breathe in and out normally through the nose, then after a relaxed exhale, pause for three to five comfortable seconds before breathing in again. Never strain — this should feel easy. It nudges your breathing rhythm back toward calm, nasal, unhurried.

WHY BOTHER?

Nasal breathing filters, warms and humidifies every breath before it reaches your lungs — everything mouth-breathing skips. The calmer your lining gets, the easier these drills become; the more you do them, the more your nose stays in practice. A rare virtuous circle.

Simplifying your routine

If you currently rely on a decongestant spray most nights, read this chapter carefully — and please loop in your doctor or pharmacist before changing anything. This is education, not medical advice.

Stopping a nightly spray abruptly usually means several very blocked nights — which is exactly why most people who try cold turkey give up. The approach many clinicians suggest instead is gradual:

- **One nostril at a time.** Keep using the spray in one nostril only, letting the other recover first. When the free side breathes comfortably on its own, switch.
- **Stretch the intervals.** If it's every night, aim for every second night on the recovering side — supported by your reset session, the saline rinse and the humidifier.
- **Expect a bumpy fortnight.** Rebound swelling fades on a scale of days to a couple of weeks. It's temporary, and it's the mechanism letting go — not your nose getting worse.

NON-NEGOTIABLE

Talk to your doctor or pharmacist first — especially if you've used a spray daily for months or years. They see this constantly (it has a name, after all) and can tailor the taper to you.

Travel, care & troubleshooting

Traveling with it

The device is pocket-sized and flight-friendly. Hotel rooms and cabins have brutally dry air — travel is when you'll be most glad you kept the routine. Pack the charging cable; a full charge typically covers a normal trip.

Weekly care, thirty seconds

- Wipe the silicone tips with warm soapy water or an alcohol wipe; dry fully.
- Store it dry, on the nightstand — visible means remembered.
- Charge when the indicator asks; don't leave it plugged in for days.

Troubleshooting

- **"I feel nothing at all."** Normal — this isn't a sensation product. Judge mornings at night 14, not warmth at night 2.
- **"The tips feel uncomfortable."** They're sitting too deep. Rest them just inside the nostril opening.
- **"I keep forgetting sessions."** Your anchor habit is too vague. Pick the exact moment ("after I brush my teeth") and leave the device next to your toothbrush tonight.
- **"No change by night 14."** Check: full 3 minutes? Nightly, no double-skips? Bedroom humidity? Still flat — write to us; the 90-day guarantee exists for a reason.

Safety notes & when to see a professional

Red light at this wavelength and intensity is non-UV, doesn't heat tissue, and has been used in photobiomodulation research for decades. Still — a few common-sense rules.

- Use it as directed: tips just inside the nostrils, ~3 minutes, once or twice daily. More is not better.
- Don't stare directly into the LEDs at close range, and don't use the device on broken or bleeding skin.
- Pregnant, nursing, managing a diagnosed condition, or under a doctor's care for anything nasal or respiratory? Check with your doctor before starting — thirty seconds of caution, always worth it.
- Keep it away from small children; the tips are small parts.

See a professional promptly if...

- The blockage is **always on the same side**, or worsening over weeks.
- There's recurring bleeding, pain, or fully lost sense of smell.
- Congestion came with an injury, or with symptoms that feel systemic (fever, severe headache).

None of those are "push through it" situations — they deserve real medical eyes, and this guide will still be here afterwards.

Three minutes tonight

Everything in this guide compresses into one sentence: **do your three minutes tonight, anchor it to a habit you already have, and judge the trend at night 14.**

No willpower project. No five-step ritual. No little bottle running your evening. Just a small red glow at the end of the day, quietly asking nothing of you.

We read every email that comes in — the glowing ones and the honest ones. If you're stuck, if something's unclear, or if night 14 arrives and you're not convinced, write to us. A real person answers, and the 90-day guarantee means the risk stays on our side of the table.

Here's to quiet nights and mornings that start clear.
— The Kalmely Team

SUPPORT

hello@kalmely.com · Replies within one business day.

Disclaimer

The Kalmely device is a wellness device intended for general well-being, relaxation and nighttime comfort. It is not a medical device and is not intended to diagnose, treat, cure or prevent any disease or medical condition. The information in this guide is provided for educational purposes only and is not medical advice; it is not a substitute for consultation with a qualified healthcare professional. Always consult your doctor or pharmacist before making changes to any medication or treatment routine, including nasal sprays, and before starting any new wellness routine if you are pregnant, nursing, or managing a diagnosed condition.

References to published research (including Neuman & Finkelstein, 1997, *Annals of Allergy, Asthma & Immunology*) describe studies of 660nm intranasal light in general and do not constitute a claim that this specific device was tested in those studies, nor a guarantee of individual results. Individual experiences vary.

© 2026 Kalmely. All rights reserved. This guide is a complimentary gift for Kalmely customers and is not for resale or redistribution.